


Expression Elite Set

Precision, Comfort, Better Outcomes



Magnification Loupe (without power bank)	1 unit
iExpressor forceps VMF-100	3 pieces
iFlo Dry Eye Mask VEM-200	2 pieces
Disposable sheet for Eye Mask	20 pieces

Soften Before You Express




By gently warming the eyelids with the iFlo Eye Mask, the oil is liquefied and the surrounding eyelid tissue becomes more relaxed, making the glands easier and safer to express. This allows for more complete clearance of blocked glands with less pressure, reducing the risk of tissue trauma.

The package includes disposable sterilized sheets, allowing iFlo Eye Mask to be hygienically shared between different patients.

Enhanced Visualization for Precise Expression

With integrated illumination and high-quality 2.5X magnification, the loupe provides a bright, enlarged view that ensures subtle blockages and gland structures are easily visible. This allows for more precise and controlled expression.




Targeted Expression from the Gland Base

The iExpressor Forceps are engineered to enable clinicians to access the base of the meibomian glands for more effective expression. The elongated, wider paddles provide secure contact with the eyelid, minimizing slippage and reducing patient discomfort during the procedure.

Success in Action: Dr. Sran Shares Her Experience

Dr. Sran shares how the **expression set** has revolutionized her workflow. With these tools, optometrists can perform **MGD expression immediately after RF or IPL treatment**—all without moving the patient from the treatment chair to the slit lamp.




VisuScience Meditech Co., Ltd.

Add: Room 501,503 No. 99 Furongzhongsan Rd.
Wuxi, Jiangsu, China

E-mail: info@visuScience.com

Web: www.visuScience.com

Distributor



Scan for Details